

兔抗 PRKCD(Phospho-Tyr64) 多克隆抗体

中文名称：兔抗 PRKCD(Phospho-Tyr64) 多克隆抗体

英文名称： Anti-PRKCD(Phospho-Tyr64) rabbit polyclonal antibody

别名： MAY1; PKCD; CVID9; nPKC-delta

抗原： PRKCD(Phospho-Tyr64)

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play distinct roles in cells.
Applications:	WB
Name of antibody:	PRKCD(Phospho-Tyr64)

Immunogen:	Peptide sequence around phosphorylation site of Tyrosine 64 (H-I-Y(p)-E-G) derived from Human PKCD.
Full name:	protein kinase C, delta
Synonyms :	MAY1; PKCD; CVID9; nPKC-delta
SwissProt:	Q05655
WB Predicted band size:	78 kDa
WB Positive control:	HepG2 cells lysates treated with PMA
WB Recommended dilution:	500-1000

