

兔抗 PRKCD 多克隆抗体

中文名称：兔抗 PRKCD 多克隆抗体

英文名称： Anti-PRKCD rabbit polyclonal antibody

别名： MAY1; PKCD; nPKC-delta

相关类别： 一抗

抗原： PRKCD

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play distinct roles in cells. The protein encoded by this gene is one of the PKC family members. Studies both in human and mice demo

	nstrate that this kinase is involved in B cell signaling and in the regulation of growth, apoptosis, and differentiation of a variety of cell types. Alternatively spliced transcript variants encoding the same protein have been observed.
Applications:	ELISA, WB, IHC
Name of antibody:	PRKCD
Immunogen:	Fusion protein of human PRKCD
Full name:	protein kinase C, delta
Synonyms :	MAY1; PKCD; nPKC-delta
SwissProt:	Q05655
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	78 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	200-1000



