

兔抗 PRKCSH 多克隆抗体

中文名称：兔抗 PRKCSH 多克隆抗体

英文名称： Anti-PRKCSH rabbit polyclonal antibody

别名： protein kinase C substrate 80K-H; GIIB; PCLD; PLD1; G19P1; PCLD1; PKCSH; AGE-R2; VASAP-60

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： PRKCSH

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes the beta-subunit of glucosidase II, an N-linked glycan-processing enzyme in the endoplasmic reticulum. The encoded protein is an acidic phosphoprotein known to be a substrate for protein kinase C. Mutations in this gene have been associated with the autosomal dominant polycystic liver disease. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	PRKCSH

Immunogen:	Fusion protein of human PRKCSH
Full name:	protein kinase C substrate 80K-H
Synonyms:	GIIB; PCLD; PLD1; G19P1; PCLD1; PKCSH; AGE-R2; VASAP-60
SwissProt:	P14314
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	59 kDa
WB Positive control:	Mouse brain tissue lysate
WB Recommended dilution:	500-2000



