

## 兔抗 PRKRA 多克隆抗体

中文名称：兔抗 PRKRA 多克隆抗体

英文名称： Anti-PRKRA rabbit polyclonal antibody

别名： protein activator of interferon induced protein kinase EIF2AK2; RAX; PACT; DYT16; HSD14

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PRKRA

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein kinase activated by double-stranded RNA which mediates the effects of interferon in response to viral infection. Mutations in this gene have been associated with dystonia. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PRKRA
<b>Immunogen:</b>	Fusion protein of human PRKRA

<b>Full name:</b>	protein activator of interferon induced protein kinase EI F2AK2
<b>Synonyms:</b>	RAX; PACT; DYT16; HSD14
<b>SwissProt:</b>	O75569
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human tonsil
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	34 kDa
<b>WB Positive control:</b>	Mouse testis tissue, Hela cell, TM4 cell, HepG2 cell lysates
<b>WB Recommended dilution:</b>	500-2000



