

兔抗 PRLHR 多克隆抗体

中文名称: 兔抗 PRLHR 多克隆抗体

英文名称: Anti-PRLHR rabbit polyclonal antibody

别名: prolactin releasing hormone receptor; GR3; GPR10; PrRPR

抗原: PRLHR

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	PRLHR is a 7-transmembrane domain receptor for prolactin-releasing hormone (PRLH; MIM 602663) that is highly expressed in anterior pituitary (Ozawa et al., 2002 [PubMed 11923475]).
Applications:	ELISA, IHC
Name of antibody:	PRLHR
Immunogen:	Synthetic peptide of human PRLHR
Full name:	prolactin releasing hormone receptor
Synonyms:	GR3; GPR10; PrRPR

SwissProt:	P49683
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	25-100

