

## 兔抗 PRPH 多克隆抗体

中文名称：兔抗 PRPH 多克隆抗体

英文名称： Anti-PRPH rabbit polyclonal antibody

别名： NEF4; PRPH1

储存： 冷冻（-20℃）

抗原： PRPH

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a cytoskeletal protein found in neurons of the peripheral nervous system. The encoded protein is a type III intermediate filament protein with homology to other cytoskeletal proteins such as desmin, and is a different protein than the peripherin found in photoreceptors. Mutations in this gene have been associated with susceptibility to amyotrophic lateral sclerosis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PRPH
<b>Immunogen:</b>	Fusion protein of human PRPH

<b>Full name:</b>	peripherin
<b>Synonyms :</b>	NEF4; PRPH1
<b>SwissProt:</b>	P41219
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human colon cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

