

兔抗 PRRG3 多克隆抗体

中文名称: 兔抗 PRRG3 多克隆抗体

英文名称: Anti-PRRG3 rabbit polyclonal antibody

别名: TMG3; PRGP3

相关类别: 一抗

抗原: PRRG3

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:

PRRG3 (proline-rich gamma-carboxyglutamic acid protein 3), also known as TMG3 (transmembrane gamma-carboxy glutamic acid protein 3) or PRGP3, is a 231 amino acid single-pass type I membrane protein found in kidney, heart and brain. PRRG3 contains an N-terminal Gla (gamma-carboxy-glutamate) domain, two PPXY motifs and is encoded by a gene that maps to human chromosome Xq28. Chromosome X consists of about 153 million base pairs and nearly 1,000 genes. Color blindness, hemophilia and Duchenne muscular dystrophy are well known X chromosome-li

	nked conditions which affect males more frequently, as males carry a single X chromosome.
Applications:	ELISA, WB
Name of antibody:	PRRG3
Immunogen:	Full length fusion protein
Full name:	proline rich Gla (G-carboxyglutamic acid) 3 (transmembrane)
Synonyms :	TMG3; PRGP3
SwissProt:	Q9BZD7
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	26 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	200-1000

