

兔抗 PRSS1 多克隆抗体

中文名称: 兔抗 PRSS1 多克隆抗体

英文名称: Anti-PRSS1 rabbit polyclonal antibody

别名: TRP1; TRY1; TRY4; TRYP1

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: PRSS1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a trypsinogen, which is a member of the trypsin family of serine proteases. This enzyme is secreted by the pancreas and cleaved to its active form in the small intestine. It is active on peptide linkages involving the carboxyl group of lysine or arginine. Mutations in this gene are associated with hereditary pancreatitis. This gene and several other trypsinogen genes are localized to the T cell receptor beta locus on chromosome 7.
Applications:	WB
Name of antibody:	PRSS1

Immunogen:	Fusion protein of human PRSS1
Full name:	protease, serine, 1 (trypsin 1)
Synonyms :	TRP1; TRY1; TRY4; TRYP1
SwissProt:	P07477
WB Predicted band size:	27 kDa
WB Positive control:	Liver tissue
WB Recommended dilution:	500-2000

