

## 兔抗 PRSS38 多克隆抗体

中文名称: 兔抗 PRSS38 多克隆抗体

英文名称: Anti-PRSS38 rabbit polyclonal antibody

别名: MPN2

抗原: PRSS38

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The function of this protein remains unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PRSS38
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	protease, serine, 38
<b>Synonyms :</b>	MPN2
<b>SwissProt:</b>	A1L453
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human cervical cancer

IHC Recommend dilution:	25-100
-------------------------	--------

