

兔抗 PSKH1 多克隆抗体

中文名称：兔抗 PSKH1 多克隆抗体

英文名称： Anti-PSKH1 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

抗 原： PSKH1

宿 主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	May be a SFC-associated serine kinase (splicing factor compartment-associated serine kinase) with a role in intranuclear SR protein (non-snRNP splicing factors containing a serine/arginine-rich domain) trafficking and pre-mRNA processing.
Applications:	WB
Name of antibody:	PSKH1
Immunogen:	Synthesized peptide derived from internal of human KPSH1.
Full name:	protein serine kinase H1
Synonyms :	0
SwissProt:	P11801
WB Predicted band size:	48 kDa

WB Positive control:	A549 cells lysate
WB Recommended dilution:	500-3000

