

兔抗 PSMA4 多克隆抗体

中文名称: 兔抗 PSMA4 多克隆抗体

英文名称: Anti-PSMA4 rabbit polyclonal antibody

别名: proteasome subunit alpha 4; HC9; PSC9; HsT17706

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PSMA4

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a core alpha subunit of the 20S proteasome, which is a highly ordered ring-shaped structure composed of four rings of 28 non-identical subunits. Proteasomes cleave peptides in an ATP- and ubiquitin-dependent manner.
Applications:	ELISA, WB, IHC
Name of antibody:	PSMA4
Immunogen:	Fusion protein of human PSMA4
Full name:	proteasome subunit alpha 4
Synonyms:	HC9; PSC9; HsT17706

SwissProt:	P25789
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	29 kDa
WB Positive control:	NIH/3T3 cell, Mouse spleen tissue, K562 cell, Hela cell, 231 cell lysates
WB Recommended dilution:	500-2000



