

兔抗 PSME1 多克隆抗体

中文名称：兔抗 PSME1 多克隆抗体

英文名称： Anti-PSME1 rabbit polyclonal antibody

别名： proteasome activator subunit 1; PA28A; IFI5111; REGalpha; PA28alpha; HEL-S-129m

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： PSME1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. Proteasomes are distributed throughout eukaryot

	<p>ic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. The immunoproteasome contains an alternate regulator, referred to as the 11S regulator or PA28, that replaces the 19S regulator. Three subunits (alpha, beta and gamma) of the 11S regulator have been identified. This gene encodes the alpha subunit of the 11S regulator, one of the two 11S subunits that is induced by gamma-interferon. Three alpha and three beta subunits combine to form a heterohexameric ring. Alternative splicing results in multiple transcript variants.</p>
Applications:	ELISA, IHC
Name of antibody:	PSME1
Immunogen:	Fusion protein of human PSME1
Full name:	proteasome activator subunit 1
Synonyms:	PA28A; IFI5111; REGalpha; PA28alpha; HEL-S-129m
SwissProt:	Q06323
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-200

