

## 兔抗 KRT31 多克隆抗体

中文名称: 兔抗 KRT31 多克隆抗体

英文名称: Anti-KRT31 rabbit polyclonal antibody

别名: HA1, Ha-1, hHa1, KRTHA1

储存: 冷冻 (-20℃)

抗原: KRT31

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	KRT31
<b>Immunogen:</b>	Fusion protein of human KRT31
<b>Full name:</b>	Keratin-31

<b>Synonyms :</b>	HA1, Ha-1, hHa1, KRTHA1
<b>SwissProt:</b>	Q15323
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer and human lung cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	47 kDa
<b>WB Positive control:</b>	A549 cells
<b>WB Recommended dilution:</b>	500-2000



