

## 兔抗 KRT36 多克隆抗体

中文名称: 兔抗 KRT36 多克隆抗体

英文名称: Anti-KRT36 rabbit polyclonal antibody

别名: keratin 36; HA6; hHa6; KRTHA6

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: KRT36

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the keratin gene family. This type I hair keratin is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	KRT36
<b>Immunogen:</b>	Synthetic peptide of human KRT36
<b>Full name:</b>	keratin 36

<b>Synonyms:</b>	HA6; hHa6; KRTHA6
<b>SwissProt:</b>	O76013
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human liver cancer
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	52 kDa
<b>WB Positive control:</b>	Human skin tissue lysate
<b>WB Recommended dilution:</b>	500-2000





