

兔抗 PSPN 多克隆抗体

中文名称：兔抗 PSPN 多克隆抗体

英文名称： Anti-PSPN rabbit polyclonal antibody

别名： PSP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PSPN

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a neurotrophic factor, belonging to the GDNF family. Neurotrophic factors are important for the proper development and maintenance of the nervous system. These factors promote neuronal survival and can prevent the neuronal degeneration associated with injury, toxin exposure, or neurodegenerative disease. The encoded protein has amino acid similarity to its other family members, glial cell line-derived neurotrophic factor and neurturin. This gene product promotes the survival of ventral midbrain dopaminergic neuron

	s in culture and prevents their degeneration after 6-hydroxydopamine treatment in vivo.
Applications:	ELISA, IHC
Name of antibody:	PSPN
Immunogen:	Synthetic peptide of human PSPN
Full name:	persephin
Synonyms :	PSP
SwissProt:	O60542
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and Human ovarian cancer
IHC Recommend dilution:	25-100

