

兔抗 ZNF83 多克隆抗体

中文名称: 兔抗 ZNF83 多克隆抗体

英文名称: Anti-ZNF83 rabbit polyclonal antibody

别名: zinc finger protein 83; HPF1; ZNF816B

抗原: ZNF83

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF83
Immunogen:	Fusion protein of human ZNF83
Full name:	zinc finger protein 83
Synonyms:	HPF1; ZNF816B
SwissProt:	P51522
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human colorectal cancer
IHC Recommend dilution:	40-200



