

兔抗 PTBP1 多克隆抗体

中文名称：兔抗 PTBP1 多克隆抗体

英文名称： Anti-PTBP1 rabbit polyclonal antibody

别名： PTB; PTB2; PTB3; PTB4; pPTB; HNRPI; PTB-1; PTB-T;

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： PTBP1

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA-binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm.
Applications:	WB, IHC

Name of antibody:	PTBP1
Immunogen:	Fusion protein of human PTBP1
Full name:	polypyrimidine tract binding protein 1
Synonyms :	PTB; PTB2; PTB3; PTB4; pPTB; HNRPI; PTB-1; PTB-T; HNRNPI; HNRNP-I
SwissProt:	P26599
IHC positive control:	Human prostate cancer and Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	57 kDa
WB Positive control:	Hela, hepg2 and HT-29 cell
WB Recommended dilution:	200-1000



