

兔抗 SMAD2(Phospho-Ser250) 多克隆抗体

中文名称：兔抗 SMAD2(Phospho-Ser250) 多克隆抗体

英文名称： Anti-SMAD2(Phospho-Ser250) rabbit polyclonal antibody

别名： JV18; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SMAD2(Phospho-Ser250)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways.
Applications:	WB
Name of antibody:	SMAD2(Phospho-Ser250)
Immunogen:	Peptide sequence around phosphorylation site of serine250(E-L-S(p)-P-T) derived from Human Smad 2 .

Full name:	SMAD family member 2
Synonyms :	JV18; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2
SwissProt:	Q15796
WB Predicted band size:	52 kDa
WB Positive control:	293 cells lysates treated with PMA
WB Recommended dilution:	500-1000

