

## 兔抗 PTCD3 多克隆抗体

中文名称: 兔抗 PTCD3 多克隆抗体

英文名称: Anti-PTCD3 rabbit polyclonal antibody

别名: pentatricopeptide repeat domain 3; MRP-S39

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PTCD3

反应种属: Human


标记物: Unconjugate



克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Mitochondrial RNA-binding protein that has a role in mitochondrial translation.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PTCD3
<b>Immunogen:</b>	Fusion protein of human PTCD3
<b>Full name:</b>	pentatricopeptide repeat domain 3
<b>Synonyms:</b>	MRP-S39
<b>SwissProt:</b>	Q96EY7
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human esophagus cancer and Human tonsil
<b>IHC Recommend dilution:</b>	100-200
<b>WB Predicted band size:</b>	79 kDa
<b>WB Positive control:</b>	K562 cell, BGC-823 Cell, HepG2 cell lysates
<b>WB Recommended dilution:</b>	500-2000





**mlbio** 酶联生物  
Good elisakit producers