

兔抗 ZNF213 多克隆抗体

中文名称: 兔抗 ZNF213 多克隆抗体

英文名称: Anti-ZNF213 rabbit polyclonal antibody

别名: CR53; ZSCAN53; ZKSCAN21

抗原: ZNF213

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	C2H2 zinc finger proteins, such as ZNF213, have bipartite structures in which one domain binds DNA or RNA and the other modulates target gene expression.
Applications:	ELISA, IHC
Name of antibody:	ZNF213
Immunogen:	Fusion protein of human ZNF213
Full name:	zinc finger protein 213
Synonyms :	CR53; ZSCAN53; ZKSCAN21
SwissProt:	O14771
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human esophagus cancer

IHC Recommend dilution: 25-100



