

兔抗 ZNF217 多克隆抗体

中文名称: 兔抗 ZNF217 多克隆抗体

英文名称: Anti-ZNF217 rabbit polyclonal antibody

别名: ZABC1

抗原: ZNF217

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

ZNF217, also known as ZABC1, is a zinc finger protein belonging to the Krüppel C2H2-type zinc-finger protein family. It localizes to the nucleus and may play a role in transcriptional repression of a variety of genes through the recruitment of co-repressor complexes containing proteins such as CtBP, HDAC1 and HDAC2. In addition, ZNF217 participates in cell differentiation and appears to function as an oncogene. Expression of ZNF217 is amplified in various tumors and overexpression of the protein can attenuate apoptotic signals and lead to epithelial cell immortalization. Binds to the p

	romoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.
Applications:	ELISA, IHC
Name of antibody:	ZNF217
Immunogen:	Fusion protein of human ZNF217
Full name:	zinc finger protein 217
Synonyms :	ZABC1
SwissProt:	O75362
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer
IHC Recommend dilution:	25-100

