

兔抗 PTGR2 多克隆抗体

中文名称: 兔抗 PTGR2 多克隆抗体

英文名称: Anti-PTGR2 rabbit polyclonal antibody

别名: PGR2; ZADH1; HEL-S-298

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: PTGR2

反应种属: Human Mouse Rat

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes an enzyme involved in the metabolism of prostaglandins. The encoded protein catalyzes the NADPH-dependent conversion of 15-keto-prostaglandin E2 to 15-keto-13,14-dihydro-prostaglandin E2. This protein may also be involved in regulating activation of the peroxisome proliferator-activated receptor. Alternative splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	PTGR2
Immunogen:	Synthesized peptide derived from internal of human ZADH

	1.
Full name:	prostaglandin reductase 2
Synonyms :	PGR2; ZADH1; HEL-S-298
SwissProt:	Q8N8N7
WB Predicted band size:	38 kDa
WB Positive control:	Huvec cells, K562 cells and HT-29 cells lysates
WB Recommended dilution:	500-3000

