

兔抗 PTGS2 多克隆抗体

中文名称：兔抗 PTGS2 多克隆抗体

英文名称： Anti-PTGS2 rabbit polyclonal antibody

别名： COX2; COX-2; PHS-2; PGG/HS; PGHS-2; hCox-2; GRIPGHS

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PTGS2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Prostaglandin-endoperoxide synthase (PTGS), also known as cyclooxygenase, is the key enzyme in prostaglandin biosynthesis, and acts both as a dioxygenase and as a peroxidase. There are two isozymes of PTGS: a constitutive PTGS1 and an inducible PTGS2, which differ in their regulation of expression and tissue distribution. This gene encodes the inducible isozyme.
Applications:	ELISA, IHC
Name of antibody:	PTGS2
Immunogen:	Fusion protein of human PTGS2

Full name:	prostaglandin-endoperoxide synthase 2 (prostaglandin G/H synthase and cyclooxygenase)
Synonyms :	COX2; COX-2; PHS-2; PGG/HS; PGHS-2; hCox-2; GRIPGHS
SwissProt:	P35354
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150

