

兔抗 ZNF256 多克隆抗体

中文名称：兔抗 ZNF256 多克隆抗体

英文名称：Anti-ZNF256 rabbit polyclonal antibody

别名：zinc finger protein 256; BMZF3; BMZF-3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZNF256

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Transcriptional repressor that plays a role in cell proliferation. Requires TRIM28 for its activity.
Applications:	ELISA, IHC
Name of antibody:	ZNF256
Immunogen:	Fusion protein of human ZNF256
Full name:	zinc finger protein 256
Synonyms:	BMZF3; BMZF-3
SwissProt:	Q9Y2P7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution: 100-500

