

兔抗 PTHLH 多克隆抗体

中文名称：兔抗 PTHLH 多克隆抗体

英文名称： Anti-PTHLH rabbit polyclonal antibody

别名： HHM, PLP, BDE2, PTHR, PTHRP

抗原： PTHLH

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the parathyroid hormone family. This hormone regulates endochondral bone development and epithelial-mesenchymal interactions during the formation of the mammary glands and teeth. This hormone is involved in lactation possibly by regulating the mobilization and transfer of calcium to the milk. The receptor of this hormone, PTHR1, is responsible for most cases of humoral hypercalcemia of malignancy. Four alternatively spliced transcript variants encoding two distinct isoforms have been observed. There is also evidence

	nce for alternative translation initiation from non-AUG (CUG and GUG) start sites, in-frame and downstream of the initiator AUG codon, to give rise to nuclear forms of this hormone.
Applications:	ELISA, IHC
Name of antibody:	PTH LH
Immunogen:	Fusion protein of human PTH LH
Full name:	parathyroid hormone-like hormone
Synonyms :	HHM, PLP, BDE2, PTHR, PTHRP
SwissProt:	P12272
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

