

## 兔抗 PTH 多克隆抗体

中文名称：兔抗 PTH 多克隆抗体

英文名称： Anti-PTH rabbit polyclonal antibody

别名： PTH1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： PTH

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a hormone secreted by parathyroid cells. This hormone elevates blood Ca <sup>2+</sup> level by dissolving the salts in bone and preventing their renal excretion. Defects in this gene are a cause of familial isolated hypoparathyroidism (FIH).
<b>Applications:</b>	WB
<b>Name of antibody:</b>	PTH
<b>Immunogen:</b>	Synthesized peptide derived from human Parathyroid Hormone antibody.

<b>Full name:</b>	parathyroid hormone
<b>Synonyms :</b>	PTH1
<b>SwissProt:</b>	P01270
<b>WB Predicted band size:</b>	13 kDa
<b>WB Positive control:</b>	A549 cells lysate
<b>WB Recommended dilution:</b>	500-3000

