

兔抗 ZNF266 多克隆抗体

中文名称：兔抗 ZNF266 多克隆抗体

英文名称：Anti-ZNF266 rabbit polyclonal antibody

别名：HZF1

相关类别：一抗

抗原：ZNF266

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a protein containing many tandem zinc-finger motifs. Zinc fingers are protein or nucleic acid-binding domains, and may be involved in a variety of functions, including regulation of transcription. This gene is located in a cluster of similar genes encoding zinc finger proteins on chromosome 19. Alternative splicing results in multiple transcript variants for this gene.
Applications:	ELISA, IHC
Name of antibody:	ZNF266
Immunogen:	Fusion protein of human ZNF266
Full name:	zinc finger protein 266

Synonyms :	HZF1
SwissProt:	Q14584
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and human esophagus cancer
IHC Recommend dilution:	25-100



