

兔抗 SMARCA4 多克隆抗体

中文名称：兔抗 SMARCA4 多克隆抗体

英文名称：Anti-SMARCA4 rabbit polyclonal antibody

别名：BRG1, SNF2, SWI2, MRD16, RTPS2, BAF190, SNF2L4, SNF2LB, hSNF2b, BAF190A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SMARCA4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the SWI/SNF family of proteins and is similar to the brahma protein of *Drosophila*. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI, which is required for transcriptional activation of genes normally repressed by chromatin. In addition, this protein can bind BRCA1, as

	well as regulate the expression of the tumorigenic protein CD44. Mutations in this gene cause rhabdoid tumor predisposition syndrome type 2. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	SMARCA4
Immunogen:	Synthetic peptide of human SMARCA4
Full name:	SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member
Synonyms :	BRG1, SNF2, SWI2, MRD16, RTPS2, BAF190, SNF2L4, SNF2LB, hSNF2b, BAF190A
SwissProt:	P51532
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	185 kDa
WB Positive control:	Hela and Jurkat cells
WB Recommended dilution:	500-2000

