

## 兔抗 PTK2 (Phospho-Tyr861) 多克隆抗体

- 中文名称：兔抗 PTK2 (Phospho-Tyr861) 多克隆抗体
- 英文名称：Anti-PTK2 (Phospho-Tyr861) rabbit polyclonal antibody
- 别名：FAK; FADK; FAK1; FRNK; PPP1R71; p125FAK; pp125FAK
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：PTK2 (Phospho-Tyr861)
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Non-receptor protein-tyrosine kinase implicated in signaling pathways involved in cell motility, proliferation and apoptosis. Activated by tyrosine-phosphorylation in response to either integrin clustering induced by cell adhesion or antibody cross-linking, or via G-protein coupled receptor (GPCR) occupancy by ligands such as bombesin or lysophosphatidic acid, or via LDL receptor occupancy. Plays a potential role in oncogenic transformations resulting in increased kinase activity.
<b>Applications:</b>	WB

<b>Name of antibody:</b>	PTK2 (Phospho-Tyr861)
<b>Immunogen:</b>	Synthetic peptide of human PTK2 (Phospho-Tyr861)
<b>Full name:</b>	protein tyrosine kinase 2 (Phospho-Tyr861)
<b>Synonyms :</b>	FAK; FADK; FAK1; FRNK; PPP1R71; p125FAK; pp125FAK
<b>SwissProt:</b>	Q05397
<b>WB Predicted band size:</b>	125 kDa
<b>WB Positive control:</b>	HT29 cells treated with EGF
<b>WB Recommended dilution:</b>	500-1000

