

## 兔抗 ZNF300 多克隆抗体

中文名称：兔抗 ZNF300 多克隆抗体

英文名称：Anti-ZNF300 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：ZNF300

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a C2H2-type zinc finger DNA binding protein and likely transcriptional regulator. The function of this protein is not yet known. Three transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ZNF300
<b>Immunogen:</b>	Fusion protein of human ZNF300
<b>Full name:</b>	zinc finger protein 300
<b>SwissProt:</b>	Q96RE9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and human liver cancer
<b>IHC Recommend dilution:</b>	25-100



