

兔抗 SMC4 多克隆抗体

中文名称：兔抗 SMC4 多克隆抗体

英文名称： Anti-SMC4 rabbit polyclonal antibody

别名： CAPC; CAP-C; SMC-4; SMC4L1; hCAP-C

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SMC4

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Members of the structural maintenance of chromosomes, or SMC, family (e.g., SMC1A; MIM 300040) are critical for mitotic chromosome condensation in frogs and for DNA repair in mammals.
Applications:	ELISA, IHC
Name of antibody:	SMC4
Immunogen:	Synthetic peptide of human SMC4
Full name:	structural maintenance of chromosomes 4
Synonyms :	CAPC; CAP-C; SMC-4; SMC4L1; hCAP-C
SwissProt:	Q9NTJ3
ELISA Recommended dilution:	1000-2000

IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

