

兔抗 ZNF329 多克隆抗体

中文名称：兔抗 ZNF329 多克隆抗体

英文名称：Anti-ZNF329 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）避光

宿 主：Rabbit

抗 原：ZNF329

反应种属：Human

标 记 物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	WB
Name of antibody:	ZNF329
Immunogen:	Synthesized peptide derived from internal of human ZNF329.
Full name:	zinc finger protein 329
Synonyms :	0
SwissProt:	Q86UD4
WB Predicted band size:	62 kDa
WB Positive control:	Huvec cells lysate
WB Recommended dilution:	500-3000

