

兔抗 ZNF337 多克隆抗体

中文名称：兔抗 ZNF337 多克隆抗体

英文名称： Anti-ZNF337 rabbit polyclonal antibody

储存： 冷冻（-20℃） 避光

抗原： ZNF337

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a zinc finger domain containing protein. The function of this protein has yet to be determined. Alternative splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	ZNF337
Immunogen:	Synthesized peptide derived from internal of human ZNF337.
Full name:	zinc finger protein 337
Synonyms :	0
SwissProt:	Q9Y3M9
WB Predicted band size:	87 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000

