

## 兔抗 SMN1 多克隆抗体

中文名称：兔抗 SMN1 多克隆抗体

英文名称： Anti-SMN1 rabbit polyclonal antibody

别名： SMA; SMN; SMA1; SMA2; SMA3; SMA4; SMA@; SMNT; BCD541; GEMIN1; TDRD16A; T-BCD541

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： SMN1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is part of a 500 kb inverted duplication on chromosome 5q13. This duplicated region contains at least four genes and repetitive elements which make it prone to rearrangements and deletions. The repetitiveness and complexity of the sequence have also caused difficulty in determining the organization of this genomic region. The telomeric and centromeric copies of this gene are nearly identical and encode the same protein.
<b>Applications:</b>	ELISA, WB

<b>Name of antibody:</b>	SMN1
<b>Immunogen:</b>	Synthetic peptide of human SMN1
<b>Full name:</b>	survival of motor neuron 1, telomeric
<b>Synonyms :</b>	SMA; SMN; SMA1; SMA2; SMA3; SMA4; SMA@; SMN T; BCD541; GEMIN1; TDRD16A; T-BCD541
<b>SwissProt:</b>	Q16637
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	32 kDa
<b>WB Positive control:</b>	Hepg2, hela and 293T cells
<b>WB Recommended dilution:</b>	200-1000

