

兔抗 PTK6 多克隆抗体

中文名称：兔抗 PTK6 多克隆抗体

英文名称： Anti-PTK6 rabbit polyclonal antibody

别名： BRK

相关类别： 一抗

抗原： PTK6

储存： 冷冻 (-20℃)

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a cytoplasmic nonreceptor protein kinase which may function as an intracellular signal transducer in epithelial tissues. Overexpression of this gene in mammary epithelial cells leads to sensitization of the cells to epidermal growth factor and results in a partially transformed phenotype. Expression of this gene has been detected at low levels in some breast tumors but not in normal breast tissue. The encoded protein has been shown to undergo autophosphorylation. Alternative splicing results in multiple transcript variants.

Applications:	ELISA, IHC
Name of antibody:	PTK6
Immunogen:	Fusion protein of human PTK6
Full name:	protein tyrosine kinase 6
Synonyms:	BRK
SwissProt:	Q13882
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	50-200

