

兔抗 ZNF365 多克隆抗体

中文名称: 兔抗 ZNF365 多克隆抗体

英文名称: Anti-ZNF365 rabbit polyclonal antibody

别名: UAN; Su48; ZNF365D

储存: 冷冻 (-20℃)

抗原: ZNF365

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes several isoforms which have different expression patterns and functions. Mutation in this gene is associated with uric acid nephrolithiasis (UAN). Alternatively spliced variants, encoding distinct proteins, have been identified.
Applications:	ELISA, WB
Name of antibody:	ZNF365
Immunogen:	Synthetic peptide of human ZNF365
Full name:	zinc finger protein 365
Synonyms :	UAN; Su48; ZNF365D
SwissProt:	Q70YC4

ELISA Recommended dilution:	1000-2000
WB Predicted band size:	24 kDa
WB Positive control:	Human liver cancer tissue
WB Recommended dilution:	200-1000



