

## 兔抗 PTK7 多克隆抗体

中文名称: 兔抗 PTK7 多克隆抗体

英文名称: Anti-PTK7 rabbit polyclonal antibody

别名: protein tyrosine kinase 7 (inactive); CCK4; CCK-4

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PTK7

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the receptor protein tyrosine kinase family of proteins that transduce extracellular signals across the cell membrane. The encoded protein lacks detectable catalytic tyrosine kinase activity, is involved in the Wnt signaling pathway and plays a role in multiple cellular processes including polarity and adhesion. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PTK7

<b>Immunogen:</b>	Synthetic peptide of human PTK7
<b>Full name:</b>	protein tyrosine kinase 7 (inactive)
<b>Synonyms:</b>	CCK4; CCK-4
<b>SwissProt:</b>	Q13308
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

