

## 兔抗 PTOV1 多克隆抗体

中文名称：兔抗 PTOV1 多克隆抗体

英文名称：Anti-PTOV1 rabbit polyclonal antibody

别名：ACID2; PTOV-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PTOV1

反应种属：Human, Mouse, Rat


标记物：Unconjugate



克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a protein that was found to be overexpressed in prostate adenocarcinomas. The encoded protein was found to interact with the lipid raft protein flotillin-1 and shuttle it from the cytoplasm to the nucleus in a cell cycle dependent manner. Alternative splicing of this gene results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PTOV1
<b>Immunogen:</b>	Synthetic peptide of human PTOV1
<b>Full name:</b>	prostate tumor overexpressed 1

<b>Synonyms :</b>	ACID2; PTOV-1
<b>SwissProt:</b>	Q86YD1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	47 kDa
<b>WB Positive control:</b>	Rat brain tissue and Mouse brain tissue
<b>WB Recommended dilution:</b>	200-1000





**mlbio** 酶联生物  
Good elisakit producers