

兔抗 ZNF410 多克隆抗体

中文名称：兔抗 ZNF410 多克隆抗体

英文名称：Anti-ZNF410 rabbit polyclonal antibody

别名：APA1; APA-1

相关类别：一抗

抗原：ZNF410

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	ZNF410 is the transcription factor that activates transcription of matrix-remodeling genes such as MMP1 during fibroblast senescence.
Applications:	ELISA, IHC
Name of antibody:	ZNF410
Immunogen:	Fusion protein of human ZNF410
Full name:	zinc finger protein 410
Synonyms:	APA1; APA-1
SwissProt:	Q86VK4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human thyroid cancer
IHC Recommend dilution:	30-150



