

## 兔抗 ZNF446 多克隆抗体

中文名称: 兔抗 ZNF446 多克隆抗体

英文名称: Anti-ZNF446 rabbit polyclonal antibody

别名: ZSCAN30; ZSCAN52; ZKSCAN20

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: ZNF446

反应种属: Human

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	May be involved in transcriptional regulation.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ZNF446
<b>Immunogen:</b>	Synthesized peptide derived from internal of human ZNF446.
<b>Full name:</b>	zinc finger protein 446
<b>Synonyms :</b>	ZSCAN30; ZSCAN52; ZKSCAN20
<b>SwissProt:</b>	Q9NWS9
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	Huvec cells lysate
<b>WB Recommended dilution:</b>	500-3000

