

兔抗 ZNF449 多克隆抗体

中文名称：兔抗 ZNF449 多克隆抗体

英文名称：Anti-ZNF449 rabbit polyclonal antibody

别名：ZSCAN19

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZNF449

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear protein that likely functions as a transcription factor. The protein includes a n N-terminal SCAN domain, and seven C2H2-type zinc finger motifs.
Applications:	ELISA, IHC
Name of antibody:	ZNF449
Immunogen:	Full length fusion protein
Full name:	zinc finger protein 449
Synonyms:	ZSCAN19
SwissProt:	Q6P9G9
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer and human breast cancer
IHC Recommend dilution:	25-100



