

## 兔抗 LAG3 多克隆抗体

中文名称：兔抗 LAG3 多克隆抗体

英文名称：Anti-LAG3 rabbit polyclonal antibody

别名：CD223

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LAG3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and chromosomal localization all indicate a close relationship of LAG4 to
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	LAG3
<b>Immunogen:</b>	Synthetic peptide of humanLAG3
<b>Full name:</b>	lymphocyte-activation gene 3
<b>Synonyms :</b>	CD223

<b>SwissProt:</b>	P18627
<b>ELISA Recommended dilution</b> :	2000-5000
<b>WB Predicted band size:</b>	57 kDa
<b>WB Positive control:</b>	Human ovarian cancer tissue
<b>WB Recommended dilution:</b>	500-2000

