

兔抗 SMYD3 多克隆抗体

中文名称：兔抗 SMYD3 多克隆抗体

英文名称： Anti-SMYD3 rabbit ployclonal antibody

别名： KMT3E; ZMYND1; ZNFN3A1; bA74P14.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SMYD3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	SMYD3
Immunogen:	Fusion protein of human SMYD3
Full name:	SET and MYND domain containing 3
Synonyms :	KMT3E; ZMYND1; ZNFN3A1; bA74P14.1
SwissProt:	Q9H7B4

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

