

兔抗 PTPN1(Phospho-Ser50) 多克隆抗体

中文名称：兔抗 PTPN1(Phospho-Ser50) 多克隆抗体

英文名称： Anti-PTPN1(Phospho-Ser50) rabbit polyclonal antibody

别名： PTP1B

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： PTPN1(Phospho-Ser50)

反应种属： Human Mouse Rat Monkey

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is the founding member of the protein tyrosine phosphatase (PTP) family, which was isolated and identified based on its enzymatic activity and amino acid sequence. PTPs catalyze the hydrolysis of the phosphate monoesters specifically on tyrosine residues. Members of the PTP family share a highly conserved catalytic motif, which is essential for the catalytic activity.
Applications:	WB
Name of antibody:	PTPN1(Phospho-Ser50)

Immunogen:	Peptide sequence around phosphorylation site of Serine 50 (D-V-S(p)-P-F) derived from Human PTP1B.
Full name:	protein tyrosine phosphatase, non-receptor type 1
Synonyms :	PTP1B
SwissProt:	P18031
WB Predicted band size:	50 kDa
WB Positive control:	COS7 cells lysates untreated or treated with UV
WB Recommended dilution:	500-1000

