

### 兔抗 ZNF467 多克隆抗体

中文名称: 兔抗 ZNF467 多克隆抗体

英文名称: Anti-ZNF467 rabbit polyclonal antibody

别名: EZI; Zfp467

储存: 冷冻 (-20℃)

抗原: ZNF467

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Transcription factor that promotes adipocyte differentiation and suppresses osteoblast differentiation in the bone marrow. Enhances the osteoclast-supporting ability of stromal cells. Binds with STAT3 the consensus sequence 5'-CTTCTGGGAAGA-3' of the acute phase response element (APRE). Transactivates several promoters including FOS, OSM and PPARG. Recruits a histone deacetylase complex (By similarity).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ZNF467
<b>Immunogen:</b>	Synthetic peptide of human ZNF467
<b>Full name:</b>	zinc finger protein 467

<b>Synonyms:</b>	EZI; Zfp467
<b>SwissProt:</b>	Q7Z7K2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and human breast cancer
<b>IHC Recommend dilution:</b>	25-100



