

兔抗 TRIB1 多克隆抗体

中文名称: 兔抗 TRIB1 多克隆抗体

英文名称: Anti-TRIB1 rabbit polyclonal antibody

别名: C8FW; GIG2; TRB1; GIG-2; SKIP1; TRB-1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: TRIB1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Interacts with MAPK kinases and regulates activation of MAP kinases. May not display kinase activity.
Applications:	ELISA, IHC
Name of antibody:	TRIB1
Immunogen:	Synthetic peptide of humanTRIB1
Full name:	tribbles pseudokinase 1
Synonyms :	C8FW; GIG2; TRB1; GIG-2; SKIP1; TRB-1
SwissProt:	Q96RU8
ELISA Recommended dilution:	1000-2000

IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	10-50

