

兔抗 ZNF473 多克隆抗体

中文名称：兔抗 ZNF473 多克隆抗体

英文名称：Anti-ZNF473 rabbit polyclonal antibody

别名：ZN473; HZFP100

相关类别：一抗

抗原：ZNF473

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Involved in histone 3'-end pre-mRNA processing by associating with U7 snRNP and interacting with SLBP /pre-mRNA complex. Increases histone 3'-end pre-mRNA processing but has no effect on U7 snRNP levels, when overexpressed. Required for cell cycle progression from G1 to S phases.
Applications:	ELISA, IHC
Name of antibody:	ZNF473
Immunogen:	Fusion protein of human ZNF473
Full name:	zinc finger protein 473
Synonyms :	ZN473; HZFP100

SwissProt:	Q8WTR7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	25-100



