

兔抗 ZNF486 多克隆抗体

中文名称：兔抗 ZNF486 多克隆抗体

英文名称： Anti-ZNF486 rabbit polyclonal antibody

别名： KRBO2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ZNF486

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. The majority of zinc-finger proteins contain a Krüppel-type DNA binding domain and a KRAB domain, which is thought to interact with KAP1, thereby recruiting histone modifying proteins. As a member of the Krüppel C2H2-type zinc-finger protein family, ZNF486 (Zinc finger protein 486), also known as KRAB domain only protein 2, is a 216 amino acid nuclear protein that contains one KRAB domain and two C2H2-type zinc fingers. The gene encoding ZNF486 maps to human chromosome 19, which consists of over 63 million bases, housing

	es approximately 1,400 genes and is recognized for having the greatest gene density of the human chromosomes. May be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF486
Immunogen:	Fusion protein of human ZNF486
Full name:	zinc finger protein 486
Synonyms :	KRBO2
SwissProt:	Q96H40
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and human thyroid cancer
IHC Recommend dilution:	25-100



