

兔抗 TRIB3 多克隆抗体

中文名称：兔抗 TRIB3 多克隆抗体

英文名称： Anti-TRIB3 rabbit polyclonal antibody

别名： NIPK; SINK; TRB3; SKIP3; C20orf97

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TRIB3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a putative protein kinase that is induced by the transcription factor NF-kappaB. The encoded protein is a negative regulator of NF-kappaB and can also sensitize cells to TNF- and TRAIL-induced apoptosis. In addition, this protein can negatively regulate the cell survival serine-threonine kinase AKT1. Differential promoter usage and alternate splicing result in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	TRIB3

Immunogen:	Full length fusion protein
Full name:	tribbles pseudokinase 3
Synonyms :	NIPK; SINK; TRB3; SKIP3; C20orf97
SwissProt:	Q96RU7
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human prostate cancer
IHC Recommend dilution:	25-100

