

兔抗 TRIM13 多克隆抗体

中文名称：兔抗 TRIM13 多克隆抗体

英文名称： Anti-TRIM13 rabbit polyclonal antibody

别名： CAR; LEU5; RFP2; DLEU5; RNF77

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： TRIM13

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This gene is located on chromosome 13 within the minimal deletion region for B-cell chronic lymphocytic leukemia. Multiple alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TRIM13

Immunogen:	Fusion protein of human TRIM13
Full name:	tripartite motif containing 13
Synonyms :	CAR; LEU5; RFP2; DLEU5; RNF77
SwissProt:	O60858
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human ovarian cancer and human brain
IHC Recommend dilution:	25-100

